

LIFETIME WARRANTY

Mega Motion offers a lifetime warranty on the steel lift mechanism. Mega Motion warrants that the steel lift mechanism on your lift chair will be free from defects under normal use and service for as long as the chair is in use by the original owner. This lifetime warranty covers the steel lift mechanism against metal fatigue, defective welds, and bushings.

NOTE: Mega Motion reserves the right to replace only the part of the steel lift mechanism that may be defective.

ONE-YEAR LIMITED WARRANTY

(FC201, LC100, LC200, LC400, LC500, SC-10, SC-11, WG325, WG701, and WG801)

For one (1) year from the date of purchase, Mega Motion will repair or replace at our option to the original purchaser any of the following parts found upon examination by an authorized representative of Mega Motion to be defective in material and/or workmanship:

- Motor
- Hand Control
- Scissor Mechanisms
- Transformer/Control Boxes

ONE-YEAR LIMITED WARRANTY

For one (1) year from the date of purchase, Mega Motion will provide labor, service calls, shipping, and/or any other charges incurred for repair of the product. All transportation costs and shipping damage incurred while submitting parts for repair or replacement are the responsibility of Mega Motion. Please contact your authorized provider for information on service visits.

This warranty covers service necessitated by product failure during normal usage. The warranty specifically does not cover failure due to customer misuse. The manufacturer reserves the right to determine if the cause of the failure is due to misuse. Customer abuse/misuse will void the warranty.

Service is not extended to fabric, buttons, or wear and tear items or issues. Labor does not cover product maintenance and/or adjustments. Only manufacturing defects are covered.

RECONDITIONED UNITS WARRANTY

All reconditioned units are covered by a six-month warranty from Mega Motion effective from the date of purchase.

WARRANTY EXCLUSIONS

This warranty does not extend to fabric and excludes steel frame chairs.

- Circumstances beyond the control of Mega Motion
- Repairs and/or modifications made to any part without specific consent from Mega Motion

Exclusions also include components with damage caused by:

- Contamination
- Abuse, misuse, accident, or negligence
- Battery fluid spillage or leakage
- Commercial use, or use other than normal
- Improper operation, maintenance, or storage



IMPLIED WARRANTIES

Implied warranties, including those of merchantability and fitness for a particular purpose, are limited to one (1) year from the date of purchase and to the extent permitted by law. Any and all other implied warranties are excluded. This is the exclusive remedy. Liabilities for consequential damages under any and all warranties are excluded.

Some states do not allow limitations on how long an implied warranty lasts or do not allow the exclusion or limitation of incidental or consequential damages. The above limitation or exclusion may not apply to you.

This warranty gives you specific rights, and you may also have other rights which vary from state to state.



Mega Motion, Inc.
182 Susquehanna Avenue
Exeter, PA 18643
1-800-MEGA-508
service@mega-motion.com

NOTE: If you ever lose or misplace your copy of this manual, contact us and we will be glad to send you a new one immediately.



PRODUCT INSERT

LIFT CHAIR WEIGHT CAPACITIES

Your lift chair was developed for your convenience. See the below specification table for information on weight capacity. Refer to the owner's manual for any additional information concerning your lift chair.

MODEL	WEIGHT CAPACITY
LC100, LC200, FC201, WG801, SC-11	325 lbs. (147 kg)
LC400, WG325, WG701, SC-10	375 lbs. (170 kg)
LC500	500 lbs. (227 kg)



MANDATORY! Stay within the specified weight capacity of your lift chair. Exceeding the weight capacity voids your lift chair warranty.



PRODUCT INSERT

LIFT CHAIR WALL MOUNT

EXTERNAL TRANSFORMER OPTIONAL WALL MOUNT DIAGRAM

Hold the transformer on the wall where you intend to mount it, and either mark the holes or drill (as shown) to secure the transformer to the wall.

